

CapsuleGuard[®]

Irrigation/Aspiration Handpiece

IA SINGLE USE

Excellent for Use in All Phases of I/A

AVAILABLE IN **3** TIP STYLES



SINGLE-USE
INSTRUMENT

STERILE 



PRODUCT FEATURES

- Flexible Tip
- Two Irrigation Ports
- One Piece Silicone Construction
- Flexible Irrigation Tube
- Soft Aspiration Port
- Flexible Sleeve Conforms to Incision
- Enhanced Visualization from Semi-transparent Silicone
- No Sharp Edges for Ease of Insertion Into the Eye

CapsuleGuard[®] IA

LEARN MORE AT bauschsurgical.com



Available in three tip styles and adaptable within a large range of incision sizes.

CapsuleGuard[®] I/A Stellaris[®] & Millennium[®]

Features a flexible tip design for cortical removal. Recommended for use in a sub 2.0mm to 2.8mm incision, the smooth irrigation and aspiration ports eliminate sharp edges for reduced risk of capsule rupture. The silicone tip design facilitates cortex removal, capsule polishing, viscoelastic removal, and IOL manipulation in the capsular bag, and the semitransparent silicone sleeve provides enhanced visualization.



Designed for use with the Stellaris[®] system

MODEL NUMBER	TIP	SIZE
85910ST	45°	Standard incision 2.2–2.8mm
85915ST	45°	MICS™, sub 2.0mm incision
85913ST	Curved	Standard incision 2.2–2.8mm
85912ST	Straight	Standard incision 2.2–2.8mm
85914ST	Straight	MICS™, sub 2.0mm incision



Designed for use with the Millennium[®] and most competitive systems

MODEL NUMBER	TIP	SIZE
85910S	45°	Standard incision 2.2–2.8mm
85915S	45°	MICS™, sub 2.0mm incision
85913S	Curved	Standard incision 2.2–2.8mm
85912S	Straight	Standard incision 2.2–2.8mm
85914S	Straight	MICS™, sub 2.0mm incision

STORZ[®]
Ophthalmic Instruments
by Bausch + Lomb

www.StorzEye.com

INSTRUMENT EXPERTS

Our Surgical Instrument Specialists are at your service to support all of your instrumentation requirements. Bausch + Lomb Storz[®] Instruments set the standard in Ophthalmic instrumentation for innovative design, quality and service. Our instruments are available directly through Bausch + Lomb Sales and Service professionals.

Our Surgical Instrument Specialists, Instrument DirectConnect Specialists and our Customer Service team are trained to Partner with your Practice to ensure you have access to the best products and the best support possible.

Customer Service (800) 338-2020

Instrument Sales Specialist (800) 325-9505