

SCLERAL FIXATION FORCEPS

 pinnacle³⁶⁰[™]

Curved shaft enables anterior access from the pars plana

Scleral Fixation Forceps

25G 11.46.25PIN =5  STERILE 

- Serrated grasping platform engages sutures without shredding
- 25ga cannula compatible with suture engaged
- Glare reducing finish for better visualization



Engineered for

 pinnacle³⁶⁰[™]

Please see [bauschsurgical.com](https://www.bauschsurgical.com) for Important Safety Information.

BAUSCH+LOMB



Synergetics[®]