


5th Edition

BAUSCH + LOMB

Single-Use Catalog

Select Pack
Custom Packs
Per Procedure Trays
Knives
Cannulas
I/A Handpieces
Cataract Instruments
Vitreoretinal Instruments
Bipolar
Trophines
Eye Shields
Other Accessories



For United States Only



Table of Contents

Select Pack	2
Custom Packs	3
Per Procedure Trays	4
Knives	5
Cannulas	9
I/A Handpieces	13
Cataract Instruments	16
Vitreoretinal Instruments	18
Bipolar	20
Trephines	21
Eye Shields	22
Other Accessories	22

Select Pack

Simplify Your Procedure Packs with the Select Pack System

Select Pack is a complete kit of components designed to streamline and optimize your cataract or vitreoretinal procedure. Combine this kit of the most popular procedural components with your choice of CapsuleGuard® I/A handpieces, blades and cannulas to create a complete procedure kit.



Select Pack Cataract



STERILE



Each Select Pack will come with standard components plus your choice of drape. Choose the drape you want in your Select Pack from the list below:

FULL BODY LENGTH

Ophthalmic Surgical Drape
2 Fluid Pouches
Breather Bar
SMS Non-Woven Material
72 x 100 inches
16 X 16" Incise Facial Drape

Select Pack Catalog #

Incise

Oval Aperature
LYN023BSBLB

MID-SIZE

Embossed Plastic Incise
Surgical Drape

LYN030INSPA



ASC Select Pack



STERILE



A concise, essential compilation of procedural components, including the most popular ophthalmic drape, priced specifically for the modern ASC.

MID-SIZE

Embossed Plastic Incise
Surgical Drape

Select Pack Catalog #

LYN030INSPA

Embossed Plastic Aperature
Surgical Drape

LYN030APSPA

Select Pack Vit



STERILE



VITREORETINAL

Vitreoretinal Pack

Select Pack Catalog #

LYN023VTBLD

80 - Cotton Tip Wood Applicators, 3"
3 - Steri-Strip Reinforced Wound Closures
10 - Eye Spears
3 - Tube Holders with Velcro
1 - Instrument Pouch, 8 Pocket
1 - Rubberband
13 - Gauze Sponges, 4" X 4"
1 - Tie Twist
1 - Instrument Wipe

2 - White Paper Towels
1 - Mayo Stand Cover
1 - Drape, Incise, 51" X 51"
1 - Drape, EENT Split, 76" X 124"
2 - Syringe 3ml, & 25g Needle
1 - Syringe, 1ml
1 - Outer Wrap/Table Cover
2 - 30g Needles



Build Your Custom Pack with these Simple Steps

STERILE



1. Determine your component list, including Storz® ophthalmics single-use instrumentation, CapsuleGuard® I/A handpieces, cannulas and Laseredge® Plus disposable knives.
2. Contact Bausch + Lomb Customer Service to initiate a quote: 800-338-2020 or SelectPackQuotes@bausch.com.
3. A Bausch + Lomb Custom Packs Coordinator will review the request and work with you to ensure you are getting the appropriate products at the best possible procedural price.
4. Once your custom-designed pack is finalized, begin receiving your inventory and spend less on handling and inventory control issues encountered when dealing with multiple suppliers.



Sterile Drapes

10

STERILE



20

STERILE



50

STERILE



10

STERILE

FULL BODY

Ophthalmic Surgical Drape
2 fluid pouches
Breather Bar
SMS Non-Woven Material
71 x 100 inches

Incise
DRA851

Oval Aperture
DRA850A

3/4 LENGTH

Ophthalmic Surgical Drape
2 fluid pouches
Breather Bar
SMS Non-Woven Material
55 x 71 inches

DRA831

DRA830A

MID-SIZE

Ophthalmic Surgical Drape
1 fluid pouch
SMS Non-Woven Material
48 x 55 inches

DRA801

DRA800A

Ophthalmic Surgical Drape
1 fluid pouch
Embossed Plastic Material
48 x 51 inches

DRA1061

DRA1059A

EQUIPMENT COVERS

Mayo Stand Cover
Clear Plastic
24 x 40 inches

DRA7352

Screen Cover
Clear Plastic

DP5300

Catalog #

Per Procedure Trays

Single-Use Instrument Per Procedure Trays

Designed for safety, consistency and convenience. Contains high precision instrumentation for the intended procedure.
Designed to optimize OR efficiency.

Cataract



STERILE



MICS™ Cataract

Catalog #	Description
PPT049	Contains 1 each: • SUH01 - Barraquer Lid Speculum • 85915S - CapsuleGuard® I/A Handpiece MICS™ Millennium®/ Other Systems 45° • SUE02-12 - Bonn Forceps • SUE04-12 - McPherson Tying Forceps • SUE06-12 - Capsulorhexis Forceps MICS • SUF01 - Koch Stop & Chop Manipulator • SUF02 - Sinskey Hook


LaserEdge® Plus Knives



Guarded Knives for Safety



E7590 Reusable Blade Tip Handle



Clear Corneal Knives



STERILE




2.75mm, Angled, Double Bevel
2.85mm, Angled, Double Bevel
2.9mm, Angled, Double Bevel
1.3mm, Angled, Double Bevel
1.6mm, Angled, Double Bevel

GUARDED HANDLE

2.85mm, Angled, Double Bevel

Catalog #	Dimensions (mm/°)			
	A	B	C	D
LP7527	2.75	3.00	15	42°
LP7528	2.85	3.20	15	42°
LP7587A	2.90	3.10	15	42°
LP7594	1.30		14	42°
LP7595	1.60	1.70	15	42°
LP7528G	2.85	3.20	11	42°



Crescent Knives



STERILE




Straight, Single Bevel, Up
Angled, Single Bevel, Up
Angled, Single Bevel, Down
Tip Only, Angled, Single Bevel Up

GUARDED HANDLE

Angled, Single Bevel, Up

Catalog #	Dimensions (mm/°)			
	A	B	C	D
LP7500	2.00	4.30	14	Straight
LP7510	2.00	4.00	14	58°
LP7510D	2.00	4.30	14	58°
LP7510AT	2.00	4.30	14	58°
LP7510G	2.00	4.30	11	45°



Implant Knives



STERILE



4.0mm, Angled, Single Bevel, Up
5.0mm, Angled, Single Bevel, Up
3.5mm, Angled, Single Bevel, Up
3.2mm, Angled, Single Bevel, Up

Catalog #	Dimensions (mm/°)		
	A	B	C
LP7560A	4.00	45°	18
LP7561A	5.00	45°	18
LP7562	3.50	45°	18
MRS32	3.20	45°	

Knives




LRI Knives

500 Microns, Straight
550 Microns, Straight
600 Microns, Straight

Catalog #

LP8500
LP8550
LP8600



MVR Knives



STERILE



19 Gauge
20 Gauge
23 Gauge
25 Gauge


GUARDED HANDLE

20 Gauge
23 Gauge
25 Gauge

Catalog # Dimensions (mm/°)
A B

LP7519	1.38	23
LP7520	1.20	23
LP7523	0.80	23
LP7525	0.65	23

LP7520G	1.20	12
LP7523G	0.80	12
LP7525G	0.65	12



Sideport/Arrow Knives



STERILE



1.1mm, Angled
0.8mm, Angled


GUARDED HANDLE

1.1mm, Angled
0.8mm, Angled

Catalog # Dimensions (mm/°)
A B C

LP7592	1.10		23
LP7593	0.80		23

LP7592G	1.10		23
LP7593G	0.80		23



Slit Knives



STERILE



2.2mm, Angled, Double Bevel
2.4mm, Angled, Double Bevel
2.5mm, Angled, Single Bevel, Up
2.65mm, Angled, Single Bevel, Up
2.65mm, Angled, Double Bevel
2.8mm, Angled, Double Bevel
2.85mm, Angled, Single Bevel, Up
3.0mm, Angled, Single Bevel, Up
3.0mm, Angled, Double Bevel
3.2mm, Angled, Single Bevel, Up
2.75mm, Angled, Single Bevel, Up

GUARDED HANDLE


2.2mm, Angled, Single Bevel, Up
2.4mm, Angled, Single Bevel, Up
2.5mm, Angled, Single Bevel, Up
2.65mm, Angled, Double Bevel
2.65mm, Angled, Single Bevel, Up
3.0mm, Angled, Single Bevel, Up
2.85mm, Angled, Single Bevel, Up

THIN BLADE

2.2mm, Angled, Single Bevel Up
2.4mm, Angled, Single Bevel Up
2.5mm, Angled, Single Bevel Up
2.65mm, Angled, Single Bevel Up
2.75mm, Angled, Single Bevel Up

Catalog #	Dimensions (mm/°)			
	A	B	C	D
LP7548ADB	2.20	2.30	14	42°
LP7549ADB	2.40	2.50	14.5	42°
LP7550A	2.50	3.30	14	42°
LP7551A	2.65	3.70	15	42°
LP7551ADB	2.65	2.80	14.5	42°
LP7552ADB	2.80	3.20	15	42°
LP7553A	2.85	3.50	14.5	42°
LP7556A	3.00	3.60	14	42°
LP7556ADB	3.00	3.10	14.5	42°
LP7557A	3.20	3.20	15	42°
LP7559A	2.75	3.40	15	42°
LP7548AG	2.20	2.20	11	42°
LP7549AG	2.40	2.40	11	42°
LP7550AG	2.50	3.30	11	42°
LP7551ADBG	2.65	2.80	11	42°
LP7551AG	2.65	3.70	11	42°
LP7556AG	3.00	3.60	11	42°
LP7553AG	2.85	3.50	12	42°
LP7747A	2.20		14	42°
LP7749A	2.40		14	42°
LP7750A	2.50		14.5	42°
LP7751A	2.65	3.70	14.5	42°
LP7759A	2.75		14.5	42°

Knives



Stab Knives

15°
22.5°
30°
45°
Tip Only, 15°
Tip Only, 30°



GUARDED HANDLE

15°

LIMITED DEPTH

15°, 1.5mm Depth, Tip Only
15°, 3.0mm Depth, Tip Only
15°, 5.0mm Depth, Tip Only

Catalog #	Dimensions (mm/°)		
	A	B	C
LP7515	15°	1.60	11
LP7522	22.5°	1.60	11
LP7530	30°	1.60	11
LP7545	45°	1.80	11
LP7515T	15°	1.60	11
LP7530T	30°	1.60	11
LP7515G	15°	1.60	11
LP7511	15°		1.50
LP7513	15°		3.00
LP7505	15°		5.00



Trapezoidal Knives



STERILE



2.75mm x 3.2mm,
Angled
2.85mm x 3.2mm,
Angled
1.5mm x 2.0mm,
Angled
1.5mm x 1.7mm,
Angled
1.6mm x 1.8mm,
Angled
2.0mm x 2.2mm,
Angled
1.8mm x 2.2mm x
2.8mm, Angled
1.8mm x 2.2mm,
Angled

GUARDED HANDLE

2.0mm x 2.2mm,
Angled

Catalog #	Dimensions (mm/°)						
	A	B	C	D	E	F	G
LP7572	2.75	3.20	14	42°			
LP7582	2.85	3.20	15	42°			
LP7596	1.50	2.00	15	42°	2.00		
LP7599	1.50	1.99	14	40°	2.20	4.50	
LP7600	1.60	1.80	14.5	42°			
LP7601	2.00	2.20	14	42°	2.30	3.10	2.25
LP7602	1.80-2.20-2.80		14.5	42°	2.00		
LP7603	1.80	2.20	14.5	42°	2.00		
LP7601G	2.00	2.20	11	42°	2.00	3.10	2.25

Anesthesia Needles



10

STERILE



23 Gauge, Atkinson Style Bevel
Retrobulbar Needle

23 Gauge, Atkinson Style Bevel
Peribulbar Needle

25 Gauge, Atkinson Style Bevel
Retrobulbar Needle

25 Gauge, Atkinson Style Bevel
Peribulbar Needle

Catalog #

Dimensions (mm)

Length (excluding hub)

M0638

38.0

M0643

32.0

M0637

38.0

M0642

32.0

Anterior Chamber Needles



10

STERILE



20 Gauge, Angled Tip

23 Gauge, Angled Tip

25 Gauge, Angled Tip

27 Gauge, Angled Tip

30 Gauge, Angled Tip

Catalog #

Tip (mm)

Dimensions (mm)

Length (excluding hub)

M2273B

7.5

22.0

M2273G

5.5

22.0

M2639

7.5

22.0

M2273A

4.5

22.0

M2273

5.5

22.0

Cannulas

Cystotomes - Irrigating



10

STERILE



CAPSULORHEXIS

25 Gauge, Formed Angled Tip

Catalog #	Dimensions (mm)
	Length (excluding hub)

M3000G

16.0

FORMED/OTHER

27 Gauge, Curved

M3001

16.0

25 Gauge, Formed

M3610

16.0

27 Gauge, Formed

M3610B

16.0

25 Gauge, Formed Reverse

M3629

16.0

25 Gauge, Short Formed

M3627

16.0

Hydrodissection



10

STERILE



⊗



25 Gauge

Catalog #	Tip (mm)	Dimension (mm) Length (excluding hub)
-----------	----------	--

M3900

7

22.0

27 Gauge

M3911

11

22.0

27 Gauge

M3902

9

22.0

25 Gauge, J Shaped

M3906A

1.75

22.0

19 Gauge, Bishop
Harmon Irrigating

M9280

8

25.0

27 Gauge, Akahoshi

M3919

2

22.0

27 Gauge, Chang

M3904

19.0

Lacrimal



10

STERILE

⊗

25 Gauge, Curved

Catalog #	Dimensions (mm) Length (excluding hub)
-----------	---

M9277

29.0

Cannulas

Miscellaneous

		Catalog #	Dimensions (mm)
			Length (excluding hub)
 	  	27 Gauge Silicone Capsule Polisher; Silicone cap to provide squeegee effect	M8605 22.0
	 	23 Gauge Cortex Extractor	M6576 22.0

Viscoelastic

		Catalog #	Tip (mm)	Dimensions (mm)
				Length (excluding hub)
 	  	27 Gauge, Hub Facilitates Flow of Viscoelastic Material	M3927 8.0	22.0

CapsuleGuard® I/A Handpieces



Benefits of CapsuleGuard® I/A Handpieces



STERILE




- This instrument features a flexible tip design for effective cortical removal.
- The smooth irrigation and aspiration ports eliminate sharp edges for reduced risk of capsule rupture.
- The silicone tip design facilitates cortex removal, capsule polishing, viscoelastic removal and IOL manipulation in the capsular bag and the semi-transparent silicone sleeve provides visualization.
- A fully assembled single-use design is provided for consistency and convenience.

CapsuleGuard® Stellaris® System



Catalog #	Tip	Size
85912ST	Straight	Standard Incision 2.2-2.8mm
85914ST	Straight	MICS™, sub 2.0mm incision
85913ST	Curved	Standard Incision 2.2-2.8mm
85910ST	45°	Standard Incision 2.2-2.8mm
85915ST	45°	MICS™, sub 2.0mm incision


CapsuleGuard® Millennium® & Other Systems



Catalog #	Tip	Size
85912S	Straight	Standard Incision 2.2-2.8mm
85914S	Straight	MICS™, sub 2.0mm incision
85913S	Curved	Standard Incision 2.2-2.8mm
85910S	45°	Standard Incision 2.2-2.8mm
85915S	45°	MICS™, sub 2.0mm incision

I/A Handpieces

Stellaris® MICS™ I/A Handpieces



12

STERILE



MICS™ 1.8mm 45° I/A
Handpiece with Sleeve



Catalog #	Description
85795ST	1.8 mm coaxial MICS™ 0.3mm aspiration port with silicone irrigation sleeve.

Aspiration Handpieces All Systems



12

STERILE



21 Gauge Aspiration
Handpiece

22 Gauge Aspiration
Handpiece

23 Gauge Aspiration
Handpiece



Catalog #	Description
85780S	Roughened Tip. 0.8mm Outer Diameter. Curved Shaft.
85790S	Roughened Tip. 0.7mm Outer Diameter. Curved Shaft.
85792S	Roughened Tip. 0.6mm Outer Diameter. Curved Shaft.

Irrigating Choppers



12

STERILE



Nagahara Irrigating Chopper for
Biaxial MICS™



Catalog #	Description
85903ST19	19 gauge. 1.8mm biaxial MICS. Central irrigation port. Outer diameter 1.1 mm. Sharpened inner side with atraumatic tip.

Irrigation Handpieces Stellaris® System



	Catalog #	Description
21 Gauge Irrigation Handpiece	85787ST	Smooth Tip. 0.8mm Outer Diameter. Curved Shaft.
22 Gauge Irrigation Handpiece	85806ST	Smooth Tip. 0.7mm Outer Diameter. Curved Shaft.
23 Gauge Irrigation Handpiece	85807ST	Smooth Tip. 0.6mm Outer Diameter. Curved Shaft.
MICS™ Irrigating Handpiece	85902ST19	19 gauge. 1.8mm biaxial MICS™. Central irrigation port. 1.1mm Outer Diameter. Sharpened inner side with atraumatic tip.

Irrigation Handpieces Millennium® & Other Systems



	Catalog #	Description
21 Gauge Irrigation Handpiece	85787S	Smooth Tip. 0.8mm Outer Diameter. Curved Shaft.
22 Gauge Irrigation Handpiece	85806S	Smooth Tip. 0.7mm Outer Diameter. Curved Shaft.
23 Gauge Irrigation Handpiece	85807S	Smooth Tip. 0.6mm Outer Diameter. Curved Shaft.
MICS™ Irrigating Handpiece	85902S19	19 gauge. 1.8mm biaxial MICS™. Central irrigation port. 1.1mm Outer diameter. Sharpened inner side with atraumatic tip

Cataract Instruments

Cataract Forceps



	Catalog #	Description
Capsulorhexis Forceps	SUE01-12	Straight jaws. Fine sharp tips.
Capsulorhexis Forceps	SUE07-12	Shafts marked at 2.75/3.0mm and 5.5/6.0mm
MICS™ Capsulorhexis Forceps	SUE06-12	Fine angled tips for 1.8mm procedures.
MICS™ Curved Capsulorhexis Forceps	SUE08-12	Fine angled tips with a curved shaft for 1.8mm procedures
Bonn Forceps	SUE02-12	1x2 teeth 0.12mm, 3.5mm tying platform. For tissue manipulation and bulbus fixation.
Tying	SUE03-12	3.5mm tying platform. May be used for IOL manipulation and cartridge loading.
McPherson Tying Forceps	SUE04-12	Angled jaws, 4.0mm tying platform. May be used for IOL cartridge loading.
Pierse Tying Forceps	SUE05-12	0.25mm notched straight jaws. 4.0mm tying platform. For tissue manipulation and bulbus fixation.

Cataract Instruments

Cataract Manipulators



STERILE



Koch Stop & Chop Manipulator

Catalog #

SUF01

Description

Angled shaft with 1.5mm tip. Designed to perform the "Stop and Chop" technique.

Sinskey Hook

SUF02

Angled shaft 10mm. Tip length 0.5mm. Diameter 0.2mm. For lens dialing and manipulation.

Phaco Spatula

SUF03

11mm spatula length. For nucleus manipulation.

Drysdale Nucleus Manipulator

SUF04

Flattened paddle tip with generous surface for manipulating tissue, rotating and cracking the nucleus.

Horizontal Manipulator with Y-Hook

SUF05

Y-shaped Nucleus Rotator with horizontal Y-hook rotator tip.

Vertical Manipulator with Y-Hook

SUF06

Y-shaped Nucleus Rotator with vertical Y-hook rotator tip.

Specula



Barraquer Lid Speculum

Catalog #

SUH01

Description

Wire speculum with closed valves. 14mm blade length.


Eye Speculum

SUH04

Adjustable Eye Speculum 16mm solid blade length.

Vitreoretinal Instruments

Forceps



12

STERILE



Catalog #	Tip	Gauge	Hub Color
SUI03-12	Asymmetric Peeling	23	Green
SUI05-12	Asymmetric Peeling	25	Blue
65023S-12	Eckardt Endgripping	23	Green
65025S-12	Eckardt Endgripping	25	Blue
SUI07-12	ILM	23	Green
SUI08-12	ILM	25	Blue
SUI11-12	Serrated	23	Green
SUI12-12	Serrated	25	Blue

Vitreoretinal Instruments

Backflush/Extrusion Handpieces



STERILE



Catalog #	Tip	Gauge	Ring Color
SUJ11	Blunt	23	Green
SUJ14	Blunt	25	Blue
SUJ13	Soft	23	Green
SUJ12	Soft Brush	23	Green

Specula



STERILE



	Catalog #	Description
Barraquer Lid Speculum	SUH01	Wire speculum with closed valves. 14mm blade length.
Eye Speculum	SUH04	Adjustable Eye Speculum 16mm solid blade length.

Bipolar

2 Pin Connection



18 GAUGE

18 Gauge Straight

Catalog #

E7918

Description

18 Gauge Curved

E7918 C

20 GAUGE

20/23 Gauge Tapered

E7923

23 GAUGE

23 Gauge Tapered

E7931

Designed for retinal endodiathermy.

25 GAUGE

25 Gauge Tapered

E7924

Designed for retinal endodiathermy.

25 Gauge Blunt

E7925

Designed for retinal endodiathermy.

CORD

Disposable Bipolar Cord

S2070 C

10 per box.

Lemo Connection



23 GAUGE

23 Gauge Tapered

Catalog #

E7930

Description

Designed for retinal endodiathermy.

25 GAUGE

25 Gauge Bipolar Diathermy Tapered Tip Lemo

E7920

Lightweight design offers precise hemostasis.

CORD

Stellaris® Bipolar Cord, Lemo Connector, Reusable

CX9430

1 per box. Reusable.

Trephines



STERILE



FOR IOWA P.K. PRESS PUNCH

Catalog #	Dimensions (mm)
	Diameter
Tanne Universal Disposable Trephine Blade	E3050 6.00 6.00
Tanne Universal Disposable Trephine Blade	E3050 6.25 6.25
Tanne Universal Disposable Trephine Blade	E3050 6.50 6.50
Tanne Universal Disposable Trephine Blade	E3050 6.75 6.75
Tanne Universal Disposable Trephine Blade	E3050 7.00 7.00
Tanne Universal Disposable Trephine Blade	E3050 7.25 7.25
Tanne Universal Disposable Trephine Blade	E3050 7.50 7.50
Tanne Universal Disposable Trephine Blade	E3050 7.75 7.75
Tanne Universal Disposable Trephine Blade	E3050 8.00 8.00
Tanne Universal Disposable Trephine Blade	E3050 8.25 8.25
Tanne Universal Disposable Trephine Blade	E3050 8.50 8.50
Tanne Universal Disposable Trephine Blade	E3050 8.75 8.75
Tanne Universal Disposable Trephine Blade	E3050 9.00 9.00
Tanne Universal Disposable Trephine Blade	E3050 9.25 9.25
Tanne Universal Disposable Trephine Blade	E3050 9.50 9.50

FOR STORZ® UNIVERSAL TREPHINE HANDLE

Trephine Blade	E3096 6.00 6.00
Trephine Blade	E3096 6.50 6.50
Trephine Blade	E3096 6.75 6.75
Trephine Blade	E3096 7.00 7.00
Trephine Blade	E3096 7.25 7.25
Trephine Blade	E3096 7.50 7.50
Trephine Blade	E3096 7.75 7.75
Trephine Blade	E3096 8.00 8.00
Trephine Blade	E3096 8.25 8.25
Trephine Blade	E3096 8.50 8.50
Trephine Blade	E3096 8.75 8.75
Trephine Blade	E3096 9.00 9.00
Trephine Blade	E3096 9.25 9.25
Trephine Blade	E3096 9.50 9.50

FOR FREEHAND USE

Trephine Blade Long	E3096L 7.00 7.00
Trephine Blade Long	E3096L 7.25 7.25
Trephine Blade Long	E3096L 7.50 7.50
Trephine Blade Long	E3096L 7.75 7.75
Trephine Blade Long	E3096L 8.00 8.00
Trephine Blade Long	E3096L 8.25 8.25
Trephine Blade Long	E3096L 8.50 8.50
Trephine Blade Long	E3096L 8.75 8.75
Trephine Blade Long	E3096L 9.00 9.00

Eye Shields, Corneal Shields & Other

Eye Shields



12



Universal Eye Shield, Adult

Catalog #	Dimensions (mm)	Description
-----------	-----------------	-------------

E5689 A	60 x 75	Clear plastic. Accommodates either eye.
---------	---------	---

Universal Eye Shield, Child

E5689 B	53.8 x 69	Clear plastic. Accommodates either eye.
---------	-----------	---



50



Combination Eye Garter Shield Cover and Fox Shield, Assorted Colors

E5691 A50	60.3 x 76.2	Assorted colors. Includes (1) E5692 D aluminum Fox Eye Shield and (1) Eye Garter Shield Cover. Latex Free.
-----------	-------------	--

Combination Eye Garter Shield Cover and Fox Shield, White Only

E5691 W50	60.3 x 76.2	White only. Includes (1) E5692 D aluminum Fox Eye Shield and (1) Eye Garter Shield Cover. Latex Free.
-----------	-------------	---

Corneal Shields



10

STERILE



Crouch Corneal Protector, Adult

Catalog #	Dimensions (mm)	Description
-----------	-----------------	-------------

E5699	13.5 x 25.8	For oculoplastic procedures for corneal and scleral protection. Dimensions enable protector to rest on the sclera and not the cornea. Easily removed with E1500 Bishop-Harmon Forceps.
-------	-------------	--

Pediatric Corneal Shield

E5699 P	17 x 17	For pediatric oculoplastic procedures for corneal and scleral protection. Dimensions enable protector to rest on the sclera and not the cornea. Easily removed with E1500 Bishop-Harmon Forceps.
---------	---------	--

Other Accessories



50

STERILE



Fluid Filter

Catalog #	Description
-----------	-------------

E4429 22	Flow rate 45cc/minute. Pore size 0.2 micron. Bi-directional membrane support rated at 5.3kg/cm ² (75 psi) via prefilled syringe and aspiration from vials.
----------	---

2 Pin Connection	20	DRA1061.....	3
65023S-12.....	18	DRA7352.....	3
65025S-12.....	18	DRA800A.....	3
85780S.....	14	DRA801.....	3
85787S.....	15	DRA830A.....	3
85787ST.....	15	DRA831.....	3
85790S.....	14	DRA850A.....	3
85792S.....	14	DRA851.....	3
85795ST.....	14	Drysdale Nucleus Manipulator.....	17
85806S.....	15	E3050 6.00.....	21
85806ST.....	15	E3050 6.25.....	21
85807S.....	15	E3050 6.50.....	21
85807ST.....	15	E3050 6.75.....	21
85902S19.....	15	E3050 7.00.....	21
85902ST19.....	15	E3050 7.25.....	21
85903ST19.....	14	E3050 7.50.....	21
85910S.....	13	E3050 7.75.....	21
85910ST.....	13	E3050 8.00.....	21
85912S.....	13	E3050 8.25.....	21
85912ST.....	13	E3050 8.50.....	21
85913S.....	13	E3050 8.75.....	21
85913ST.....	13	E3050 9.00.....	21
85914S.....	13	E3050 9.25.....	21
85914ST.....	13	E3050 9.50.....	21
85915S.....	13	E3096 6.00.....	21
85915ST.....	13	E3096 6.50.....	21
ASC Select Pack.....	2	E3096 6.75.....	21
Backflush/Extrusion.....	19	E3096 7.00.....	21
Barraquer Lid Speculum.....	17, 19	E3096 7.25.....	21
BIPOLAR.....	20	E3096 7.50.....	21
Bonn Forceps.....	16	E3096 7.75.....	21
CANNULAS.....	9-12	E3096 8.00.....	21
Cannulas, Anesthesia Needles.....	9	E3096 8.25.....	21
Cannulas, Anterior Chamber Needles.....	9	E3096 8.50.....	21
Cannulas, Cystotomes - Irrigating.....	10	E3096 8.75.....	21
Cannulas, Hydrodissection.....	11	E3096 9.00.....	21
Cannulas, Lacrimal.....	11	E3096 9.25.....	21
Cannulas, Miscellaneous.....	12	E3096 9.50.....	21
Cannulas, Viscoelastic.....	12	E3096L 7.00.....	21
Cannulas, Vitreoretinal.....	12	E3096L 7.25.....	21
Capsulorhexis Forceps.....	16	E3096L 7.50.....	21
Cataract Forceps.....	16	E3096L 7.75.....	21
CATARACT INSTRUMENTS.....	16, 17	E3096L 8.00.....	21
Cataract Manipulators.....	17	E3096L 8.25.....	21
Combination Eye Garter.....	22	E3096L 8.50.....	21
Cord.....	21	E3096L 8.75.....	21
Corneal Shield, Pediatric.....	22	E3096L 9.00.....	21
CORNEAL SHIELDS & OTHER.....	22	E4429 22.....	22
Crouch Corneal Protector.....	22	E5689 A.....	22
CUSTOM PACKS.....	3	E5689 B.....	22
CX9430.....	20	E5691 A50.....	22
DP5300.....	3	E5691 W50.....	22
DRA1059A.....	3	E5699.....	22

Index

E5699 P	22	LP7520G	6
E7918	20	LP7522	8
E7918 C	20	LP7523	6
E7920	20	LP7523G	6
E7923	20	LP7525	6
E7924	20	LP7525G	6
E7925	20	LP7527	5
E7930	20	LP7528	5
E7931	20	LP7528G	5
EQUIPMENT COVERS	3	LP7530	8
EYE SHIELDS	22	LP7530T	8
Eye Speculum	17, 19	LP7545	8
Fluid Filter	22	LP7548ADB	7
Forceps, Cataract	16	LP7548AG	7
Forceps, Asymmetrical Peeling	18	LP7549ADB	7
Forceps, Endgripping	18	LP7549AG	7
Forceps, ILM	18	LP7550A	7
Forceps, Serrated	18	LP7550AG	7
Handpieces, Aspiration	14	LP7551A	7
Handpieces, Backflush/Extrusion	19	LP7551ADB	7
Handpieces, CapsuleGuard® MILLENNIUM®	13	LP7551ADBG	7
Handpieces, CapsuleGuard® Stellaris® System	13	LP7551AG	7
Handpieces, Irrigation Choppers	14	LP7552ADB	7
Handpieces, Irrigation Stellaris® System	15	LP7553A	7
Handpieces, Stellaris® MICS™	14	LP7553AG	7
Horizontal Manipulator with Y-Hook	17	LP7556A	7
I/A HANDPIECES	13-15	LP7556ADB	7
KNIVES	5	LP7556AG	7
Knives, Clear Corneal	5	LP7557A	7
Knives, Crescent	5	LP7559A	7
Knives, Implant	5	LP7560A	5
Knives, LASEREDGE® Plus	5-8	LP7561A	5
Knives, LRI	6	LP7562	5
Knives, MVR	6	LP7572	8
Knives, Sideport/Arrow	6	LP7582	8
Knives, Slit	7	LP7587A	5
Knives, Stab	8	LP7592	6
Knives, Trapezoidal	8	LP7592G	6
Koch Stop & Chop Manipulator	17	LP7593	6
Lemo Connection	20	LP7593G	6
LP7500	5	LP7594	5
LP7505	8	LP7595	5
LP7510	5	LP7596	8
LP7510AT	5	LP7599	8
LP7510D	5	LP7600	8
LP7510G	5	LP7601	8
LP7511	8	LP7601G	8
LP7513	8	LP7602	8
LP7515	8	LP7603	8
LP7515G	8	LP7747A	7
LP7515T	8	LP7749A	7
LP7519	6	LP7750A	7
LP7520	6	LP7751A	7

LP7759A.....	7	Shield, Universal Adult.....	22
LP8500.....	6	Shield, Universal Child.....	22
LP8550.....	6	Sinskey Hook.....	17
LP8600.....	6	Specula.....	17
LYN023BSBLB.....	2	Sterile Drapes.....	3
LYN023VTBLD.....	2	SUE01-12.....	16
LYN030APSPA.....	2	SUE02-12.....	16
LYN030INSPA.....	2	SUE03-12.....	16
M0637.....	9	SUE04-12.....	16
M0638.....	9	SUE05-12.....	16
M0642.....	9	SUE06-12.....	16
M0643.....	9	SUE07-12.....	16
M2273.....	9	SUE08-12.....	16
M2273A.....	9	SUF01.....	17
M2273B.....	9	SUF02.....	17
M2273G.....	9	SUF03.....	17
M2639.....	9	SUF04.....	17
M3000G.....	10	SUF05.....	17
M3001.....	10	SUF06.....	17
M3610.....	10	SUH01.....	17, 19
M3610B.....	10	SUH04.....	17, 19
M3627.....	10	SUI03-12.....	18
M3629.....	10	SUI05-12.....	18
M3900.....	11	SUI07-12.....	18
M3902.....	11	SUI08-12.....	18
M3904.....	11	SUI11-12.....	18
M3906A.....	11	SUI12-12.....	18
M3911.....	11	SUJ11.....	19
M3919.....	11	SUJ12.....	19
M3927.....	12	SUJ13.....	19
M6576.....	12	SUJ14.....	19
M8605.....	12	Tray, MICS Cataract.....	4
M9277.....	11	TREPHINES.....	21
M9280.....	11	Trephines, for Iowa P.K. Press.....	21
MRS32.....	5	Trephines, for Storz® Universal.....	21
Mayo Stand Cover.....	3	Trephines, Freehand Use.....	21
McPherson Tying Forceps.....	16	Tying Forceps.....	16
MICS™ Capsulorhexis Forceps.....	16	Vertical Manipulator with Y-Hook.....	17
Needles, Anesthesia.....	9	VITREORETINAL INSTRUMENTS.....	18, 19
Needles, Anterior Chamber.....	9	Vitreoretinal Pack.....	2
Needles, Hydrodissection.....	11		
Other Accessories.....	22		
PER PROCEDURE TRAYS.....	4		
Phaco Spatula.....	17		
Pierse Tying Forceps.....	16		
PPT049.....	4		
S2070 C.....	20		
Screen Cover Clear Plastic.....	2		
SELECT PACK.....	2		
Select Pack Cataract.....	2		
Select Pack Vit.....	2		
Shields, Corneal.....	22		
Shield, Pediatric Corneal.....	22		



Contact your Instrument Representative at:

(800) 338-2020

www.StorzEye.com

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.
All other product/brand names are trademarks of their respective owners.
©2025 Bausch & Lomb Incorporated or its affiliates. SSU.0001.USA.25