

# UNIVERSAL VISCIOUS FLUID CONTROL PACK




## Efficiency and Compatibility in One Pack

The Universal Viscous Fluid Control Pack provides the utility of two packs in a single product, allowing for the efficient injection and extraction of viscous fluids.

The syringe adapter can accommodate both a pre-filled syringe and a 10cc syringe, allowing the surgeon to utilize either 5000 mPa-s or 1000 mPa-s silicone oil. High-flow injection cannulas are provided in both 23ga and 25ga, eliminating the need to enlarge sclerotomies during small-gauge procedures. To make these benefits available to as many surgeons as possible, the pack is designed to be compatible with most vitrectomy platforms.

### Universal Viscous Fluid Control Pack Includes:

- A 10cc syringe
- B Syringe plunger
- C Black piston
- D Syringe cap
- E Tubing set with syringe coupler
- F Female-to-female adapter for tubing set (for DORC® EVA and Alcon Constellation®)
- G 20ga cannula
- H 23ga High-Flow VFI cannula
- I 25ga High-Flow VFI cannula
- J Extraction cannula



- Use with 1000 or 5000 mPa-s silicone oil
- High-flow injection cannulas increase flow up to 209% for 23ga and 359% for 25ga\*
- High-flow extraction cannula seats over the hub of the instrument cannula
- Compatible with the Stellaris Elite®, Stellaris®PC, Alcon Constellation®, and DORC® EVA vitrectomy platforms

\*As compared to standard injection cannulas

### Ordering Information

| Product # | Description (6/Box)                  |
|-----------|--------------------------------------|
| BL7600    | Universal Viscous Fluid Control Pack |

**BAUSCH+LOMB**  **Synergetics®**